

Anti- Phospho-Fyn (Tyr530) Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A100538-01	50 µl
A100538-02	100 µl

Specificity	Anti- Phospho-Fyn (Tyr530) (human mouse rat)
Source	Rabbit Polyclonal
Application	WB IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P06241

Other Names: DKFZp762A223; FA-D2; FA4; FADC; FADC2; FAD; FAD2; FANCD; FANCD2; Fanconi anemia group D2 protein; Fanconi anemia, complementation group D2; FLJ23826; Protein FADC2; type 4 Fanconi pancytopenia

Specificity and Sensitivity

Phospho-Fyn (Tyr530) Antibody detects endogenous levels of Fyn only when phosphorylated at Tyrosine 530.

Source and Purification

A synthesized peptide derived from human Fyn around the phosphorylation site of Tyrosine 530.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:200


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of Fyn phosphorylation expression in H_2O_2 treated 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP