

Anti- Mammaglobin Polyclonal Antibody



Catalog No.	Size
A300124-01	50 µl
A300124-02	100 µl

Specificity	Anti- Mammaglobin (Human, Rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC IF
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: Q13296

Other Names: MGB1, MAMMAGLOBIN 1, secretoglobin, family 2A, member 2

Specificity and Sensitivity

Mammaglobin antibody detects endogenous levels of total mammaglobin protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human mammaglobin .

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:100
IF: 1:100~1:500 ELISA: 1:1000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions


Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images




Peptide - +

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using Mammaglobin Antibody.



Peptide - +

Western blot analysis of extracts from HepG2 cells, using Mammaglobin Antibody.



Peptide - +

Immunofluorescence analysis of HepG2 cells, using Mammaglobin Antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E032221: Dylight 488, Donkey Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP