

Anti- Mammaglobin B Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300125-01	50 µl
A300125-02	100 µl

Specificity	Anti- Mammaglobin B (Human)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: O75556

Other Names: LACRYGLOBIN, MGB2, MAMMAGLOBIN 2, MAMMAGLOBIN B, secretoglobin family 2A member 1

Specificity and Sensitivity

Mammaglobin B antibody detects endogenous levels of total Mammaglobin B protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human Mammaglobin B.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 IHC: 1:50~1:100
ELISA: 1:20000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Mammaglobin B antibody.



Western blot analysis of extracts from HepG2 cells, using Mammaglobin B antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E032221: Dylight 488, Donkey Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP