

Anti- MCM5 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300129-01	50 µl
A300129-02	100 µl

Specificity	Anti- MCM5 (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC IF
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P33992

Other Names: CDC46, DNA replication licensing factor MCM5, minichromosome maintenance deficient protein 5

Specificity and Sensitivity

MCM5 antibody detects endogenous levels of total MCM5 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human MCM5 antibody.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 ELISA: 1:20000

IHC: 1:50~1:100 IF: 1:100~1:500


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions


Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images




Peptide - +

Immunohistochemical analysis of paraffin-embedded human tonsil tissue using MCM5 antibody.



Peptide - +

Western blot analysis of extracts from HepG2 cells, using MCM5 antibody.



Peptide - +

Immunofluorescence analysis of HepG2 cells, using MCM5 antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E032221: Dylight 488, Donkey Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP