

Anti- MMP-23 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300140-01	50 µl
A300140-02	100 µl

Specificity	Anti- MMP-23 (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: Q13561

Other Names: Matrix metalloproteinase-23 [Precursor], EC 3.4.24.-, Matrix metalloproteinase 21, MMP-21
Matrix metalloproteinase 22, MMP-22, Femalysin, MIFR-1

Specificity and Sensitivity

MMP-23 antibody detects endogenous levels of total MMP-23 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human MMP-23 antibody.
The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 ELISA: 1:10000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of extracts from SK-OV3 cells, using MMP-23 antibody.

Related Products

- E030120 : HRP, Goat Anti-Rabbit IgG(H+L)
- E032221: Dylight 488, Donkey Anti-Rabbit IgG(H+L)
- E021010: Anti-GAPDH Mouse Monoclonal Antibody
- E021020: Anti-beta Actin Mouse Monoclonal Antibody
- E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP