

Anti- MMP-7 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300142-01	50 µl
A300142-02	100 µl

Specificity	Anti- MMP-7 (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P09237

Other Names: PUMP1, PUTATIVE METALLOPROTEINASE I, matrix metalloproteinase 7 (matrilysin, uterine)

Specificity and Sensitivity

MMP-7 antibody detects endogenous levels of total MMP-7 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human MMP-7.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 ELISA: 1:5000

IHC: 1:50~1:100

Storage Buffer


Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid


repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Peptide - +

Immunohistochemical analysis of paraffin-embedded human colon carcinoma tissue using MMP-7 antibody.



Western blot analysis of extracts from COS7 cells, using MMP-7 antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E032221: Dylight 488, Donkey Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP