

Anti- P300/CBP Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300159-01	50 µl
A300159-02	100 µl

Specificity	Anti- P300/CBP (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: Q09472/Q92793

Other Names: CREB-binding protein, CREBBP

E1A-associated protein p300, EP300, P300 transcription coactivator homolog

Specificity and Sensitivity

P300/CBP antibody detects endogenous levels of total P300/CBP protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human P300/CBP.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:100 ELISA: 1:10000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Peptide - +
Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using P300/CBP Antibody.



Peptide - +
Western blot analysis of extracts from COLO205 cells, using P300/CBP Antibody.

Related Products

- PW001: Super ECL Assay kit
- E030120 : HRP, Goat Anti-Rabbit IgG(H+L)
- E030220 : AP, Goat Anti-Rabbit IgG(H+L)
- E021010: Anti-GAPDH Mouse Monoclonal Antibody
- E021020: Anti-beta Actin Mouse Monoclonal Antibody
- E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP