

# Anti- Parkin Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300163-01	50 µl
A300163-02	100 µl

<b>Specificity</b>	Anti- Parkin (human mouse rat)
<b>Source</b>	Rabbit Polyclonal
<b>Application</b>	WB ELISA IHC
<b>Form</b>	Liquid, 1 mg/ml

## Product

**Swiss-Prot No.:** O60260

**Other Names:** PARK2, PRKN, PRKN2, Parkinson disease protein 2, Parkinson juvenile disease protein 2, ubiquitin E3 ligase PRKN

## Specificity and Sensitivity

Parkin antibody detects endogenous levels of total Parkin protein.

## Source and Purification

The antiserum was produced against synthesized peptide derived from human Parkin antibody.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000    IHC: 1:50~1:100    ELISA: 1:5000


## Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage Instructions


Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

## Images



Peptide                    -                    +

Immunohistochemical analysis of paraffin-embedded human brain tissue, using Parkin antibody.



Peptide                    -                    +

Western blot analysis of extracts from Jukat cells, using Parkin antibody.

## Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP