

Anti- PDGFR α Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300169-01	50 μ l
A300169-02	100 μ l

Specificity	Anti- PDGFR α (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC IF
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P16234

Other Names: Alpha platelet-derived growth factor receptor precursor, CD140a, PDGF-R-alpha, PDGFR-alpha, PGFRA, kinase PDGFR-alpha

Specificity and Sensitivity

PDGFR α antibody detects endogenous levels of total PDGFR α protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from internal of human PDGFR α .

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:100

IF: 1:100~1:500 ELISA: 1:20000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at $-20^{\circ}C$ from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at $4^{\circ}C$ for 3 months.


Images



Peptide - +
Immunohistochemical analysis of paraffin-embedded human brain tissue using PDGFR α antibody.



Peptide - +
Immunofluorescence analysis of HeLa cells, using PDGFR α antibody.



Peptide - +
Western blot analysis of extracts from HepG2 cells, using PDGFR α antibody.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP