

Anti-Prostate-specific Antigen Polyclonal Antibody



Catalog No.	Size
A300177-01	50 µl
A300177-02	100 µl

Specificity	Anti-Prostate-specific Antigen (human mouse)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P07288

Other Names: 3.4.21.77; APS, KLK2A1, PSA, gamma-seminoprotein, hK3, semenogelase, seminin, ANTIGEN, PROSTATE-SPECIFIC, P-30 antigen, PROSTATE-SPECIFIC ANTIGEN kallikrein 3, (prostate specific antigen), prostate specific antigen isoform 1 preproprotein prostate specific antigen isoform 2

Specificity and Sensitivity

Prostate-specific Antigen antibody detects endogenous levels of total prostate-specific Antigen protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human prostate-specific Antigen.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:100 ELISA: 1:20000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions


Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Peptide - +

Immunohistochemical analysis of paraffin-embedded human prostate carcinoma tissue using Prostate-specific Antigen antibody.



Western blot analysis of extracts from mouse heart, using Prostate-specific Antigen antibody.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP