

Anti- Synuclein-pan Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300198-01	50 µl
A300198-02	100 µl

Specificity	Anti- Synuclein-pan (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P37840

Other Names: NACP, Non-A beta component of AD amyloid, Non-A4 component of amyloid precursor, SYN, SYUA, alpha-synuclein, synuclein alpha

Specificity and Sensitivity

Synuclein-pan antibody detects endogenous levels of total Synuclein protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human synuclein.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:100 ELISA: 1:5000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Peptide - +
Immunohistochemical analysis of paraffin-embedded human brain tissue using Synuclein-pan antibody.



Western blot analysis of extracts from rat heart, using Synuclein-pan antibody.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP