

Anti-Thyroid Hormone Receptor β Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300206-01	50 μ l
A300206-02	100 μ l

Specificity	Anti- Thyroid Hormone Receptor β (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P10828

Other Names: C-ERBA-BETA, ERBA2, NR1A2, THB1, THR1

Specificity and Sensitivity

Thyroid Hormone Receptor β antibody detects endogenous levels of total Thyroid Hormone Receptor β protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human Thyroid Hormone Receptor β antibody. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 IHC: 1:50~1:100 ELISA: 1:5000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at $-20^{\circ}C$ from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at $4^{\circ}C$ for 3 months.

Images



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using Thyroid Hormone Receptor β antibody.



Western blot analysis of extracts from LOVO cells, using Thyroid Hormone Receptor β antibody.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP