

# Anti- SENP6 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300228-01	50 µl
A300228-02	100 µl

<b>Specificity</b>	Anti- SENP6 (human)
<b>Source</b>	Rabbit Polyclonal
<b>Application</b>	ELISA IHC
<b>Form</b>	Liquid, 1 mg/ml

## Product

**Swiss-Prot No.:** Q9GZR1

**Other Names:** EC 3.4.22.-, FKSG6, KIAA0797, Protease FKSG6, SEN6, SSP1, SUMO-1 specific protease 1, SUSP1, Sentrin-specific protease 6, Sentrin/SUMO-specific protease SENP6

## Specificity and Sensitivity

SENP6 antibody detects endogenous levels of total SENP6 protein.

## Source and Purification

The antiserum was produced against synthesized peptide derived from Internal of human SENP6. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

IHC: 1:50~1:100    ELISA: 1:10000


## Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

## Images



Peptide                      -                      +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using SENP6 antibody.

## Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP