

Anti-MAP2K3 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300531-01	50 µl
A300531-02	100 µl

Specificity	Anti-MAP2K3 (human mouse)
Source	Rabbit Polyclonal
Application	WB ELISA
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P46734

Other Names: Dual specificity mitogen-activated protein kinase kinase 3; MAP kinase kinase 3; MAP2K3; MAPK/ERK kinase 3; MAPKK 3; MAPKK3; MEK 3; MEK3; mitogen-activated protein kinase kinase 3; MKK3; MP2K3; PRKMK3

Specificity and Sensitivity

MAP2K3 antibody detects endogenous levels of total MAP2K3 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from N-terminal of human MAP2K3.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 ELISA: 1:10000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of extracts from 3T3 cells, using MAP2K3 antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP