

Anti-NCBP2 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300560-01	50 µl
A300560-02	100 µl

Specificity	Anti-NCBP2 (human mouse)
Source	Rabbit Polyclonal
Application	WB ELISA
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P52298

Other Names: 20 kDa nuclear cap binding protein, CBP20, NCBP 20 kDa subunit, NCBP-interacting protein 1, NCBP2, NIP1, Nuclear cap binding protein subunit 2

Specificity and Sensitivity

NCBP2 antibody detects endogenous levels of total NCBP2 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from internal of human NCBP2.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 ELISA: 1:5000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of extracts from COLO205 cells, using NCBP2 antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP