

Anti-GCVK_HHV6Z Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300562-01	50 µl
A300562-02	100 µl

Specificity	Anti-GCVK_HHV6Z (human)
Source	Rabbit Polyclonal
Application	WB ELISA
Form	Liquid, 1 mg/ml

Images

Product

Swiss-Prot No.: P52446

Other Names: Probable ganciclovir kinase, EC=2.7.1.-

Specificity and Sensitivity

GCVK_HHV6Z antibody detects endogenous levels of total GCVK_HHV6Z protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from internal of human GCVK_HHV6Z.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:


WB: 1:500~1:1000 ELISA: 1:10000

Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.



Western blot analysis of extracts from 3T3 cells, using GCVK_HHV6Z antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP