

Anti-SFRS3 Polyclonal Antibody



<u>Catalo No.</u>	<u>Size</u>
A300605-01	50 µl
A300605-02	100 µl

Specificity	Anti-SFRS3 (human mouse)
Source	Rabbit Polyclonal
Application	WB ELISA IF
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P84103

Other Names: Splicing factor, arginine/serine-rich 3, Pre-mRNA-splicing factor SRP20, SFRS3, SRP20

Specificity and Sensitivity

SFRS3 antibody detects endogenous levels of total SFRS3 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from internal of human SFRS3.

The antibody was affinity-purified by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 IF: 1:100~1:500 ELISA: 1:5000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of extracts from HeLa cells and Jurkat cells, using SFRS3 antibody.



Immunofluorescence analysis of MCF-7 cells, using SFRS3 antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP