

# Anti-CSTF1 Polyclonal Antibody



<u>Catalo No.</u>	<u>Size</u>
A300625-01	50 µl
A300625-02	100 µl

<b>Specificity</b>	Anti-CSTF1 (human mouse rat)
<b>Source</b>	Rabbit Polyclonal
<b>Application</b>	WB ELISA
<b>Form</b>	Liquid, 1 mg/ml

## Images

### Product

**Swiss-Prot No.:** Q05048

**Other Names:** CF-1 50 kDa subunit, CSTF 50 kDa subunit, CSTF1, Cleavage stimulation factor, 50 kDa subunit

### Specificity and Sensitivity

CSTF1 antibody detects endogenous levels of total CSTF1 protein.

### Source and Purification

The antiserum was produced against synthesized peptide derived from N-terminal of human CSTF1.

The antibody was affinity-purified by affinity-chromatography using epitope-specific immunogen.

### Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:


WB: 1:500~1:1000    ELISA: 1:20000

### Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

### Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.



Western blot analysis of extracts from 293 cells and HUVEC cells, using CSTF1 antibody.

### Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP