

# Anti-Myelodysplasia Syndrome 1 Polyclonal Antibody



<u>Catalo No.</u>	<u>Size</u>
A300657-01	50 µl
A300657-02	100 µl

<b>Specificity</b>	Anti-Myelodysplasia Syndrome 1 (human)
<b>Source</b>	Rabbit Polyclonal
<b>Application</b>	WB ELISA
<b>Form</b>	Liquid, 1 mg/ml

## Product

**Swiss-Prot No.:** Q13465

**Other Names:** Myelodysplasia syndrome-associated protein 1, Myelodysplasia syndrome 1 protein

## Specificity and Sensitivity

Myelodysplasia Syndrome 1 antibody detects endogenous levels of total Myelodysplasia Syndrome 1 protein.

## Source and Purification

The antiserum was produced against synthesized peptide derived from N-terminal of human Myelodysplasia Syndrome 1.

The antibody was affinity-purified by affinity-chromatography using epitope-specific immunogen.

## Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000      ELISA: 1:1000


## Storage Buffer

Rabbit IgG in phosphate buffered saline (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

## Images



Western blot analysis of extracts from 562 cells, using FOXE3 antibody.

## Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP