

Anti-TRA-2 α Polyclonal Antibody



<u>Catalo No.</u>	<u>Size</u>
A300666-01	50 μ l
A300666-02	100 μ l

Specificity	Anti-TRA-2 α (human mouse)
Source	Rabbit Polyclonal
Application	WB IHC ELISA IF
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: Q13595

Other Names: Transformer-2 protein homolog, TRA-2 alpha

Specificity and Sensitivity

TRA-2 α antibody detects endogenous levels of total TRA-2 α protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from internal of human TRA-2 α .

The antibody was affinity-purified by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:1000 IHC: 1:50~1:100

IF: 1:100~1:500 ELISA: 1:20000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.


Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TRA-2 α antibody.



Western blot analysis of extracts from HUVEC cells, using TRA-2 α antibody.



Immunofluorescence analysis of HeLa cells, using TRA-2 α antibody.

Related Products

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP