

Anti- BAX Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A301201-01	50 µl
A301201-02	100 µl

Specificity	Anti- BAX (human mouse rat)
Source	Rabbit Polyclonal
Application	IHC ELISA IF
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P55269

Other Names: BAX protein, cytoplasmic isoform delta, BAX

Specificity and Sensitivity

BAX antibody detects endogenous levels of total BAX protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from internal of human BAX.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

IF: 1:100~1:500 IHC: 1:50~1:100 ELISA: 1:10000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions


Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using BAX antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using BAX antibody.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP