

Anti-RAB5C Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A302225-01	50 µl
A302225-02	100 µl

Specificity	Anti- Rab5c (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA:
Form	Liquid, 1 mg/ml

Catalog Number: A302225

Concentration: 1 mg/ml

Swiss-Prot No.: P51148

Other Names: L1880; RAB5L; Ras-related protein Rab-5C

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RAB5C.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RAB5C antibody detects endogenous levels of total RAB5C protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:1000 ELISA: 1:5000

References:

Han H.J., Cytogenet. Cell Genet. 73:137-139(1996).

Albertsen H.M., Nat. Genet. 7:472-479(1994).

Clemens D.L., Infect. Immun. 68:2671-2684(2000).

