

# Anti-CBP Tag Mouse Monoclonal Antibody( Mix)



<u>Catalog No.</u>	<u>Size</u>
E022470-01	100 µl
E022470-02	500 µl
E022470-03	50 µl

<b>Specificity</b>	Anti- CBP Tag Mouse Monoclonal Antibody( Mix)
<b>Source</b>	Mouse Monoclonal
<b>Application</b>	WB IP
<b>Form</b>	Liquid, 1 mg/ml

## Product

**Other Names:** calmodulin binding peptide

## Specificity and Sensitivity

Antibody can detects CBP recombinant protein.

## Source and Purification

The antiserum was produced against synthesized CBP tag. The antibody was affinity-purified by affinity-chromatography using epitope-specific immunogen.

## Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:5000-10000 IP: 1:200


## Storage Buffer

Mouse IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

## Images



Western blot analysis of CBP recombinant protein with CBP Tag Mouse Monoclonal Antibody( Mix) diluted at 1:5,000 (line1) and 1:10,000 (line2).

## Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP