

## HRP AffiniPure Goat Anti-Rabbit IgG(H+L)



<u>Catalog No.</u>	<u>Size</u>
E030120-01	100µl
E030120-02	500µl

<b>Product Type</b>	Secondary Antibody
<b>Physical State</b>	Liquid
<b>Antibody Source</b>	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads
<b>Antibody Purity</b>	The antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins.
<b>Suggested Dilution Range</b>	1:500 - 1:5000 for immunohisto/cytochemistry 1:5000 - 1:100000 for ELISA and Western blotting with chromogenic substrates 1:10000 - 1:200000 for Western blotting with ECL substrates.
<b>Storage Buffer</b>	0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH7.6
<b>Stabilizer</b>	15 mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)
<b>Preservative</b>	None (Warning: Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.)
<b>Storage Instructions</b>	Store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Note</b>	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However, said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.