

DyLight 800 AffiniPure Rabbit Anti-Goat IgG(H+L)



<u>Catalog No.</u>	<u>Size</u>
E032830-01	100µl
E032830-02	500µl

Product Type	Secondary Antibody
Physical State	Liquid
Antibody Source	The antibody was affinity purified from antisera by immunoaffinity chromatography using purified goat IgG
Antibody Purity	The antibody reacts with whole molecule goat IgG. It also reacts with the light chains of other goat immunoglobulins.
Suggested Dilution Range	1:100- 1:800 for most applications
Fluorophore	DyLight 800-TFP ester Amax=770 nm; Emax=794 nm
Fluorophore/Protein Ratio	1:1-3:1
Storage Buffer	100 mM phosphate, 150 mM NaCl, 50% Glycerol pH 7.4
Stabilizer	1% Bovine Serum Albumin(IgG-Free)
Preservative	0.02% Sodium Azide
Storage Instructions	Store at -20°C. Avoid repeated freeze / thaw cycles.

Note The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However, said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

