

# 放大眼镜

## 放大倍率固定放大镜

|        | 2.5X 放大镜<br>(短距离) | 2.5X 放大镜<br>(长距离) | 3.5X 放大镜<br>(短距离) | 3.5X 放大镜<br>(长距离) |
|--------|-------------------|-------------------|-------------------|-------------------|
| Part # | LPM-P25S-3307     | LPM-P25L-4907     | LPM-P35S-3307     | LPM-P35L-4907     |
| 工作距离   | 330 +/- 50mm      | 490 +/- 50mm      | 330 +/- 50mm      | 490 +/- 50mm      |
| 视角     | 20-150°           | 20-150°           | 20-150°           | 20-150°           |
| 视野     | 100mm             | 110mm             | 60mm              | 65mm              |



2.5X Premium



3.5X Premium

## 放大倍率可变放大镜

|        | 2.5X-3.5X Variable |
|--------|--------------------|
| Part # | LPM-P35V-4407      |
| 工作距离   | 440mm-250mm        |
| 景深     | 130mm-50mm         |
| 视野     | 100mm-40mm         |



## 放大倍率可变放大镜 (夹式)

|        | 2.5X Clip-On  | 3.5X Clip-On  |
|--------|---------------|---------------|
| Part # | LPM-C25F-3807 | LPM-C35F-3807 |
| 工作距离   | 380 +/- 30mm  | 380 +/- 20mm  |
| 景深     | 60mm          | 40mm          |
| 视野     | 80mm          | 50mm          |



P 770.270.1394  
F 770.270.2389

865 Marathon Parkway  
Lawrenceville GA 30046

LWScientific.com